



Hyperpurex®

S

Smart series

Laboratory water purification system



Excellent Both Internally and Externally Efficient and Concise





Hyperpurex[®]

S smart series (SU/SD/SRS/SP)

Integrated Pure Water/Ultrapure Water System

Breakthrough design to highlight the aesthetics of science and technology.

S series, using the innovative automatic control system and LCD display, embedding new purification cartridges with patented structure, stable and reliable RO system with higher ion rejection rate, and DI ion-exchange cartridges with larger capacity, equipping with built-in 1.8-liter pressure water tank, more economical and cost-effective, is the economic choice for lab pure water.

System output ranges from 13 to 60 liters/h (SP series, with pure water inlet, output is up to 2 liters/minute). It can simultaneously produce ultrapure water (18.2MΩ.cm), high pure water (>17.5MΩ.cm) or pure water(<5μs/cm)^[1]. The quality of pure water fully meets or exceeds the requirements of water quality standard, specified by ASTM D1193-06, GB/T 11446.1-2013, GB/T 33087-2016, GB/T 6682-2008, CP, EP, USP, JP, CAP, CLSI, etc.


^[1] SU/SD series products can produce single RO water (ion rejection rate ≥ 98%). SRS series products can produce double RO water (<5μs/cm).

Application Area:

- HPLC、UPLC、LC-MS
- ICP-MS、ICP-AES、AAS、GC-MS
- MALDI-TOF-MS、IC、TOC analysis
- Electrochemical, spectrophotometric determination
- Preparation of microbial media and reagents
- Cell culture, PCR, IVF
- Protein purification, electrophoresis, biochemistry
- Proteomics, genomics, immunoassay
- Feed water of laboratory instruments, such as: autoclave, bottle washing machine, environmental test chamber, water bath, etc.

Smart and Concise System Design

Create Excellent Quality Both Internally and Externally



1 Easy-to-use automatic control system

- White LCD display, size up to 68*87mm, shows intuitive and easy-to-read running status.
- Real-time running status display of flushing, producing water, water full, water shortage and maintenance.
- 2 water quality sensors (RO water, DI water or UP water) can monitor water quality and alarm real-timely.
- 4 (PP/PC/RO/DI) consumables life management function, can real-timely display remaining life of consumables, automatically remind expiration replacement, and avoid water quality declines.
- Acrylic touch panel, with 3-button layout, achieving fast system setup, RO forced flushing and easy 2-way water dispensing function, brings efficient and convenient operating experience.



2 Powerful 12-inch pretreated cartridge

- 2 in 1 composite cartridge-high performance activated carbon fiber and deep folding membrane, accuracy of 5μm, eliminates particles and adsorbs organics and residual chlorine efficiently, to avoid carbon powder precipitation maximumly.



3 Rigorous double RO system^[1]

- Double RO system can remove up to 99% soluble inorganic ions, 99% soluble organics, microorganisms and particles.
- Compared with single RO system, the double RO water quality can be stable <5μs/cm (feed water conductivity <1500μs/cm), and the life of the ultrapure unit is longer.
- Equipped with DuPont RO membrane, to achieve combination of long life, stability and high ion rejection rate.
- Auto-flushing function of RO module, to effectively prevent scale and prolong the life of the membrane.
- The automatic discharge function of unqualified RO water can ensure that the RO water quality is suitable to enter the back-end module.
- Integral package of discarded RO module, easy to install and maintain.


^[1] Applicable to SRS.



4 High performance purification cartridge^[1]

- Patented cartridge structure uses full droop flow mode to prevent the stratification of resin and ensure the exchange capacity of cartridge.
- The resin filling capacity per cartridge is up to 1.36 liters, and up to 2 cartridges can be equipped every host, with a total filling capacity of 2.72 liters, achieving greater ion exchange capacity and significantly reducing the running cost.
- All DuPont resin and high purity material of column ensure absolute 18.2MΩ.cm of ultrapure water resistivity and reduce TOC precipitation.

^[1] According to different model, cartridge configuration is different. For details, refer to the product manual.



5 Double wavelength UV module^[1]

- Long-life ultraviolet lamp (185&254nm), combined with SUS316L flow shell, can reduce the value of TOC to ≤ 2ppb^[2], and can achieve efficient sterilization and inhibit bacterial growth, suitable for HPLC, UPLC, LC-MS and other precision instruments.

^[1] Applicable to ultrapure water systems equipped with UV module.

^[2] The values vary depending on the nature and concentration of contaminants in source water.



6 Ultrafiltration module^[1]

- With PES membrane and MWCO>5000D, effectively removes pyrogen/endotoxin, RNase, DNase, and produces nuclease-free, proteinase-free and bacterial-free ultrapure water, suitable for life science applications, such as cell culture/IVF.

^[1] Applicable to ultrapure water systems equipped with UF module.



8 MF terminal microfilter^[1]

- (0.45±0.2)μm double-layer PES membrane ensures microbial retention, effectively removes particles and bacteria, and meets critical application requirements.

^[1] Applicable to ultrapure or high-pure water system. For details, refer to the product manual.



8 Built-in 1.8-liter pressure water tank

- With dual functions of water storage and pressurization, FDA approved, its fully enclosed structure effectively isolates air, and prevent the touching of CO₂ and other pollutants with pure water. Up to 100 liters is optional volume.
- 60 or 120 liters pure water tank with liquid level sensor, equipped with air filter, is optional to achieve more professional pure water storage.



9 UF terminal ultrafilter^[1]

- With PES membrane and MWCO>15000D, effectively removes pyrogen/endotoxin, RNase, DNase, and produces nuclease-free, proteinase-free and bacterial-free ultrapure water, suitable for life science applications, such as cell culture/IVF.

^[1] Optional accessory for ultrapure water system only.

Combination Of Technology & Aesthetics

Creating highlights both inside and out



All injection molded housing

- New and advanced manufacturing process bring compelling customer experience.
- With geometric surfaces and simple lines, to show rich three-dimensional sense. With extraordinary imagination, to highlight the aesthetics of science and technology. Beautiful & Easy to use.

Innovative design of cartridge structure

- Patented 3-chamber design, compatible with packaging of PP/PC/RO/DI cartridge, to ensure consistency.
- Patented clamping mechanism, easier and more efficient to install and replace the cartridge.
- Patented error-proofing design, effective to avoid installation errors of different cartridges.
- 12-inch cylinder with 1.36L resin filling capacity brings more bigger ion exchange capacity and more effective filtration.
- Encrypted long serial number verification code can identify the authenticity of cartridges, record the use and replacement of cartridges, and ensure the safety of the system.



Powerful HiDis water dispenser arm (Optional)

- Color display, to monitor dispensing resistivity, water temperature, flow rate, single and cumulative water quantity.
- General, quantitative, instant - 3 water dispensing modes cycle, meeting with needs of different water dispensing mode.
- It can be fixed on the bracket in any direction of 360 degrees horizontally, making dispensing water more flexible in different directions.
- Function of circulating with the host can always ensure the quality of pure water.
- Equipped with 0.2μm MF terminal microfilter or UF terminal ultrafilter, to produce bacterial-free, nuclease-free, proteinase-free ultrapure water.
- Up to 5 sets of HiDis water dispenser arm can be connected to one host, fully covering the pure water usage range on the laboratory table.

Dispensing Mode	2021-03-08 11:35	
	18.2 MQ.cm @25°C	
	3.2 ppb	25.2 °C
General	Flux 0000 ml/min	
	Amount	500 ml



Professional PE pure water tank (Optional)

- Material: HDPE, double layer design. Anti-UV inhibitor is added to the outer layer to prevent the growth of algae inside and improve the durability of the tank. Pure PE raw material is used in inner layer to reduce material precipitation and ensure water quality safety.
- Drainage valve is installed at the cone bottom, which can empty the water tank and ensure thorough cleaning.
- Feeding from the bottom can reduce CO₂ absorption.
- The enlarged cover with seal can prevent air from entering and facilitate manual cleaning.
- Compound air filter is in the standard configuration, containing special packing and microporous membrane, to absorb CO₂ and organics, and filtrate bacteria and particles.
- UV disinfection module is optional to sterilize tank and inhibit the growth of bacteria in the tank.
- Equipped with an independent pressure sensor, independent level control module and LCD display, it can display the liquid level, storage (L) and storage percentage (%) of the water tank in real-time in the form of dynamic icons. A clear glance for storage status.

SU Smart series

Integrated Ultrapure Water System


—Ultrapure water, RO^{1st} water

With tap water inlet, using the innovative automatic control system and LCD display, embedding new purification cartridges with patented structure, stable and reliable single RO system with higher ion rejection rate, and DI ion-exchange cartridges with larger capacity, equipping with built-in 1.8-liter pressure water tank.

System output: 20, 40, 60 liters/h. It can simultaneously produce ultrapure water (18.2MΩ.cm) and single RO water. The quality of pure water fully meets or exceeds the requirements of water quality standard specified by ASTM D1193-06, GB/T 11446.1-2013, GB/T 33087-2016, GB/T 6682-2008, CP, EP, USP, JP, CAP, CLSI, etc.



Flow Diagram



SU Specifications

Name	Standard	Low TOC	Eliminating endotoxin	Synthesizing
Model	SU-20/40/60	SU-20/40/60UV	SU-20/40/60UF	SU-20/40/60UVF
Production rate ^[1]	20 series: 20 L/hour, 40 series: 40 L/hour, 60 series: 60 L/hour			
Dispensing rate ^[2]	Up to 2 liters/minute	Up to 2 liters/minute	Up to 2 liters/minute	Up to 2 liters/minute
Ultrapure water quality ^[3]				
Resistivity (25°C) ^[4]	18.2 MΩ.cm	18.2 MΩ.cm	18.2 MΩ.cm	18.2 MΩ.cm
Conductivity (25°C)	0.055 μs/cm	0.055 μs/cm	0.055 μs/cm	0.055 μs/cm
TOC ^[5]	5 ppb ^[6]	2 ppb ^[7]	5 ppb ^[6]	2 ppb ^[7]
Particles ^[8]	<1 /ml (>0.2μm)	<1 /ml (>0.2μm)	<1 /ml (>0.2μm)	<1 /ml (>0.2μm)
Bacteria ^[9]	<0.01 CFU/ml	<0.01 CFU/ml	<0.01 CFU/ml	<0.01 CFU/ml
Endotoxin ^[10]	N/A	N/A	<0.001 EU/ml	<0.001 EU/ml
RNases ^[10]	N/A	N/A	1 pg/ml	1 pg/ml
DNases ^[10]	N/A	N/A	5 pg/ml	5 pg/ml
Protease ^[10]	N/A	N/A	0.15 μg/ml	0.15 μg/ml
RO ^{1st} water quality ^[3]				
Ion rejection rate	98%-99% (with new RO module)	98%-99% (with new RO module)	98%-99% (with new RO module)	98%-99% (with new RO module)
Organic rejection rate	>99% (MW>300 Dalton)	>99% (MW>300 Dalton)	>99% (MW>300 Dalton)	>99% (MW>300 Dalton)
Particles and bacteria rejection rate	>99%	>99%	>99%	>99%
Feed water requirements				
Water source type	Tap water	Tap water	Tap water	Tap water
Pressure	1-6 bar	1-6 bar	1-6 bar	1-6 bar
Temperature	5-40 °C	5-40 °C	5-40 °C	5-40 °C
Conductivity	<2000 μs/cm	<2000 μs/cm	<2000 μs/cm	<2000 μs/cm
Total hardness (In CaCO ₃)	<300 ppm	<300 ppm	<300 ppm	<300 ppm
TOC	<2000 ppb	<2000 ppb	<2000 ppb	<2000 ppb
Free chlorine	<3 ppm	<3 ppm	<3 ppm	<3 ppm
PH	4-10	4-10	4-10	4-10
Dissolved CO ₂	<30 ppm	<30 ppm	<30 ppm	<30 ppm
Power supply	100-240V, 50/60Hz	100-240V, 50/60Hz	100-240V, 50/60Hz	100-240V, 50/60Hz
Total Power	20 series: 48W, 40 series: 72W, 60 series: 120W			
Dimension (L×W×H)	Main host: 273×555×568mm	Main host: 273×555×568mm	Main host: 273×555×568mm	Main host: 273×555×568mm
weight	About 21KG	About 21KG	About 21KG	About 21KG
Standard configuration	Main host 1 set All cartridges 1 set Built-in 1.8-liter water tank 1 set	Main host 1 set All cartridges 1 set Built-in 1.8-liter water tank 1 set	Main host 1 set All cartridges 1 set Built-in 1.8-liter water tank 1 set	Main host 1 set All cartridges 1 set Built-in 1.8-liter water tank 1 set

[1] Affected by inlet water quality, pressure, temperature and status of RO membrane
[2] Affected by the tank status and terminal filter
[3] The following values are typical and may vary depending on the nature and concentration of feed water contaminants
[4] According to USP requirements, the resistivity can be displayed as a non-temperature-compensated value
[5] Affected by the type of organics
[6] Inlet TOC<1000ppb, follow professional operating procedures and correct sampling conditions

[7] Inlet TOC<50ppb, follow professional operating procedures and correct sampling conditions
[8] Equip with terminal microfilter and follow professional operating procedures and correct sampling conditions
[9] Equip with terminal microfilter and follow professional operating procedures and correct sampling conditions
[10] Equip with terminal ultrafilter and follow professional operating procedures and correct sampling conditions